

## Specimen Test Menu

Collect, store, and transport all of these samples at room temperature (15-30° Celcius), unless otherwise specified.

### GYN Cytology



#### ThinPrep® Pap Test 33717-0

ThinPrep	33717-0
ThinPrep with Imaging	33717-0
High-Risk HPV (w/ reflex Genotyping if positive)	49896-4
Chlamydia / Gonorrhea (CT/GC)	64017-7
Trichomonas vaginalis (TV)	46154-1

### Aptima® Specimen Collection Kits

#### Unisex Collection Kit (Female & Male)

Chlamydia / Gonorrhea (CT/GC)	64017-7
Mycoplasma genitalium (MG) (female only)	69935-5
Trichomonas vaginalis (TV)	46154-1



#### Urine Collection Kit (Female & Male)

Chlamydia / Gonorrhea (CT/GC)	64017-7
Mycoplasma genitalium (MG)	69935-5
*Trichomonas vaginalis (TV)	46154-1

\*Not FDA approved, sent to reference laboratory



### Non GYN Cytology 33716-2

#### CytoLyt® Solution

Abdominal fluid	Fine-needle aspiration
Bladder washings	Misc. fluid
Breast aspirate	Pericardial fluid
Breast discharge	Pleural fluid
Bronchial brushings	Sputum
Bronchial washings	Synovial fluid
Esophageal brushings	Urine



#### Multitest Swab Collection Kit (Female ONLY)

Chlamydia / Gonorrhea (CT/GC)	64017-7
Mycoplasma genitalium (MG)	69935-5
Trichomonas vaginalis (TV)	46154-1
Bacterial vaginosis (BV)	69568-4
Vaginitis Candida / Trichomonas (CV/TV)	62460-1
Herpes simplex virus 1 & 2 (HSV)	80354-4



#### Multitest Swab Collection Kit (Male ONLY)

Herpes simplex virus 1 & 2 (HSV)	80354-4
----------------------------------	---------

#### Multitest Swab Collection Kit (Throat or Rectal)

Chlamydia / Gonorrhea (CT/GC)	64017-7
-------------------------------	---------

### Anatomic Pathology Surgical Pathology Studies 11529-5

#### 10% Neutral Buffered Formalin

Tissue biopsy



#### Michel's Solution

Tissue biopsy (Special studies)  
Immunofluorescence



#### RPMI

Tissue biopsy (Special studies)  
FLOW Cytometry  
\*Must be refrigerated (2-8° Celcius)



The correct specimen collection technique is very important in preserving the specimen. Clinical Pathology Associates provides the collection containers for your convenience.

Let us know if you have any questions by contacting:

Laboratory 325-670-6527

Marketing Representative 325-670-6513  
Jennifer Talley

Laboratory Manager 325-670-6522  
Vince Van Cleave

## Clinical HXs we may need for tests with for Cytology Specimens\*

ICD10s

Well Woman Exam	Z01.419
Apparent Infection (Other specified inflammation of vagina and vulva)	N76.89
Irregular Menses	N92.6
Abnormal Bleeding	N93.9
Cervicitis	N72
High risk heterosexual behavior	Z72.51
Pregnant	Z34.90
PostPartum	Z39.2

### Other ICD10s that could be included, if more specific

#### Chlamydia/Gonorrhea

Z20.2 - Contact with and (suspected) exposure to infections with a predominantly sexual mode of transmission

N73.9 - Pelvic Inflammatory disease

N34.1 - non-specific urethritis

N34.2 - other urethritis

#### Trichomonas

Z20.2 - Contact with and (suspected) exposure to infections with a predominantly sexual mode of transmission

N73.9 - Pelvic Inflammatory disease

N34.1 - non-specific urethritis

N34.2 - other urethritis

#### BV/CV

N76.0 - Acute vaginitis

N76.1 - Subacute and chronic vaginitis

O86.13 - Vaginitis following delivery

**N89.8 - Other specified noninflammatory disorders of vagina:**

Cyst of vagina

Hymen tag

Hymenal tag

Nontraumatic non-obstetric vaginal hematoma

Old vaginal laceration

Vaginal cyst

**Vaginal discharge**

Vaginal hematoma

Vaginal irritation

Vaginal mass

**N89.9 - Noninflammatory disorder of vagina, unspecified:**

Non inflammatory disorder of the vagina

Vaginal lesion

N77.1 - Vaginitis, vulvitis and vulvovaginitis in diseases classified elsewhere

#### Mycoplasma Genitalium

N72 -Inflammatory disease of cervix uteri; cervicitis with or w/o erosion

N34.1 - non-specific urethritis

N73.9 -Female pelvic inflammatory disease, unspecified

Z20.2 - Contact with and (suspected) exposure to infections with a predominantly sexual mode of transmission

#### HSV

A60.0 Herpesviral infection of genitalia and urogenital tract